

IHE Integration Statement			
Vendor	Product Name	Version	Date
Hitachi Medical Corporation	EUB-6500	V06-03A	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Acquisition Modality	Broad Worklist	
Patient Information Reconciliation (PIR)	Acquisition Modality	None	
Internet address for vendor's IHE information:			
Links to Standards Conformance Statements for the Implementation			
HL7			
DICOM			
Links to general information on the IHE			
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	

IS_B6500-01