

MR Imaging System

ECHELON OVAL

IHE Integration Statement

 **Hitachi Medical Corporation**
Tokyo, Japan

E1E-BM7411-1

Copyright © Hitachi Medical Corporation. 2011. All rights reserved.

| IHE Integration Statement | | | |
|---|---------------------------|---|--|
| Vendor | Product Name | Version | Date |
| Hitachi Medical Corporation | ECHELON OVAL | V4.0 | 1-Sep-2011 |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below: | | | |
| Integration Profiles Implemented | Actors Implemented | Options Implemented | |
| Scheduled Workflow (SWF) | Acquisition Modality | Broad Worklist Query Patient Based Worklist Query | |
| Patient Information Reconciliation (PIR) | Acquisition Modality | None | |
| Consistent Presentation of Images (CPI) | Acquisition Modality | None | |
| | Print Composer | None | |
| Basic Security (SEC) | Acquisition Modality | Encryption | |
| | Print Composer | | |
| | Secure Node | | |
| Audit Trail and Node Authentication (ATNA) | Acquisition Modality | Encryption | |
| | Print Composer | | |
| | Secure Node | | |
| Consistent Time (CT) | Time Client | None | |
| Key Image Note (KIN) | Acquisition Modality | None | |
| Portable Data for Imaging (PDI) | Portable Media Creator | None | |
| Internet address for vendor's IHE information: | | | |
| Links to Standards Conformance Statements for the Implementation | | | |
| HL7 | | | |
| DICOM | | | |
| Links to general information on the IHE | | | |
| In North America: www.ihe.net/ | | In Europe: www.ihe-europe.org/ | In Japan: www.ihe-j.org/ |