

MR Imaging System

Model AIRIS II / Altaire / AIRIS Elite

IHE Integration Statement

HITACHI MEDICAL CORPORATION

Tokyo, Japan

E1E-BM6800-1

Copyright © Hitachi Medical Corporation. 2007. All rights reserved.

IHE Integration Statement			
Vendor	Product Name	Version	Date
Hitachi Medical Corporation	AIRIS II/AIRIS Elite/Altaire	V5.1	9-May-2007
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Acquisition Modality	Broad Worklist Query Patient Based Worklist Query	
Patient Information Reconciliation (PIR)	Acquisition Modality	None	
Consistent Presentation of Images (CPI)	Acquisition Modality	None	
	Print Composer	None	
Basic Security (SEC)	Acquisition Modality	Encryption	
	Print Composer		
	Secure Node		
Consistent Time (CT)	Time Client	None	
Internet address for vendor's IHE information:			
Links to Standards Conformance Statements for the Implementation			
HL7			
DICOM			
Links to general information on the IHE			
In North America: www.ihe.net		In Europe: www.ihe-europe.org	In Japan: www.ihe-j.org